



VARIABLE FLOW PERISTALTIC PUMPS

Model D-25Vplus

Code 1.9749.00 and 1.9749.50



MANUAL

May 2023

Marked 

DINTER^{sa}

c/ Encarnació, 123 -125. Tel. +34 932 846 962. Fax +34 932 104 307
e-mail: dinter@dinko.es www.dinko.es 08024 - Barcelona

INDEX

	Page
1- GENERAL INTRODUCTION.....	3
2- PACKING LIST.....	3
3- RECEPTION	4
4- HEAD DESCRIPTION 253-3r.....	5
5- START UP.....	8
6- CHANGE OF TUBES.....	11
7- ORDERING INFORMATION.....	12
8- MAINTENANCE - SPARE PARTS.....	13
9- ACCESSORIES.....	13
10- CHANGING THE FUSES.....	15
11- WARRANTY	16
12- EC DECLARATION OF CONFORMITY	17
13- OTHER EQUIPMENT.....	18

1- GENERAL INTRODUCTION

Peristaltic pumps pump all kinds of liquid substances without coming into contact with mechanical elements as in other pumps.

They are easy to use with minimal maintenance.

The pumped substance is impelled into an elastic tube by the vacuum generated by rotors that successively press and release the surface of the tube.

The liquid passes directly from its container to another without any contamination and without going back when the pump stops, since the tube remains pressed by the roller.

Some aggressive substances prevent the use of conventional pumps and make peristaltic pumps very useful for the transfer or dispensing of such substances.

A large number of different size tubes or hoses made of materials resistant to various hazardous substances are available.

The following instructions are intended to ensure correct reception and use of the device, and the safety of the user.

To this end, we recommend reading this manual in detail before proceeding to unpack the device and subsequent use.

For proper conservation of the device, it is necessary to avoid its installation in areas with corrosive atmospheres or exposed to liquid splashes.

Avoid using the device when there is the possibility of generating explosive and flammable gas mixtures.

2- PACKING LIST

Description	Code	Quantity
Peristaltic Pump D-25Vplus	1. 9749.00 / 1.9749.50	1
Set of connections		1
Power cord		1
Instruction Manual		1
Warranty		1

3- RECEPTION

To ensure correct reception, use of the device, and user safety, we recommend reading this manual in detail before proceeding to unpack the device and subsequent use, and especially the following points:

3.1-THE MANUAL

This manual must be permanently kept within the equipment user's reach.

3.2-UNPACKING

Unpack the appliance, checking that the contents match the packing list.
Immediately notify any eventuality.

3.3-EXPLOSIVE MIXTURES

Avoid using the device when there is the possibility of generating explosive gas mixtures and flammable. The ATEX Directive is not covered.

3.4-LIABILITY






According to European regulation 89/655/CEE, the lack of adequate maintenance and the alteration of component, exempts the manufacturer from any responsibility for any damage that may occur.

3.5-REPAIRS

Devices to be sent to *DINKO technical services* must be **clean and disinfected** .
Otherwise, they will be rejected and returned with postage paid by the owner.

3.6-SIGNS AND SYMBOLS

Pay attention at all times to the danger warning signs and symbols that will appear in this manual or on labels attached to the body of the Pump such as those shown below.

SIGN/ SYMBOL	INTERPRETATION-MEANING
	Avoid finger contact with moving parts
	Danger-Risk-Caution
Before opening DISCONNECT the network cable Before remove cover PULL OUT plug 	Before accessing the interior of the Pump, disconnect the power cable from the network
	Possible overheating - Do not touch
110-230V AC 50/60Hz	AC power supply voltage
110V AC 60Hz	AC power supply voltage
12V DC or 24V DC	DC power supply voltage
	Disposal of waste electrical and electronic equipment by users within the European Union. It is not disposable as household waste. Deliver to the agency for recycling of electronic equipment. Contact your local office, the store where you purchased the equipment, or your household waste disposal service. Recycling helps conserve natural resources. Make sure it is recycled protecting human health and the environment.

4- HEAD DESCRIPTION 253-3r

The D-25Vplus peristaltic pumps in this manual mount the 253-3r head that allows easy access to the tube for removal when it must be replaced.

In order to change the tube, the top of the head must be removed. To be able to extract it, the lever (number 1 in the photo) must always be located to the left of the head. Once located, it must be extracted upwards with the help of the two retainers located on both sides (number 3 in the photo).

Once the upper part has been extracted, we proceed to place the tube with which we want to work, taking care that it is centered in the yellow guides located on both sides (number 8 in the photo).


Now we are going to replace the upper part, check that the lever (number 1 in the photo) is located on the left. Insert the upper part carefully and making sure that the 4 rails (number 7 in the photo) that have the removable part and the piece in the central part (number 9 in the photo) enter their corresponding guides. Once correctly oriented, press down until the two parts of the head, the fixed one (number 5 in the photo) and the removable one (number 6 in the photo) are at the same level and we check that the tube has been well centered. Once we reach this position, we move the lever (number 1 in the photo), from the left to the right, slowly, checking that all the parts are well aligned.

When removing the tube, slide the buttons (3) up to move the tube retainers. Install the new tube ensuring that it is well centered and proceed to move the lever to the right.

Lightly smearing the peristaltic tube with high-density silicone grease improves its life.

NEVER CARRY OUT THIS OPERATION WITH THE EQUIPMENT RUNNING.

4.1 HEAD 253-3r:



The 253-3r head accepts a complementary head, which is placed in front of the already installed head. In order to install it, it is necessary to remove the protector 2 to face the complementary head with the axis and insert the two assembly screws.

With 2 heads it is possible to cancel the typical cadence of the peristaltic flow. To do this, when locating the complementary head, the rotor rollers must be oriented in opposition to the main head rollers.

Each head will have a flow rate, depending on the diameter of the pipe placed. If they are the same diameter, the flow rate will double.

The option of using a Y connection to join the tubes of both heads will allow the cancellation of the peristaltic pulse. It must be remembered that the final discharge tubes and especially the suction ones must be of a larger diameter than the head tube. If this is not possible, the total flow will be somewhat less than the theoretical expected.

4.2- DESCRIPTION OF THE FRONT PANEL



A-Digital Reader

1-Key direction of rotation

2-Key direction of rotation


3-Key Full

4-Start / Stop key-Memo start

5 and 6- Decrease and increase keys.

7- Switch for programming

4.3- DESCRIPTION OF THE REAR PANEL (photo b).



1- Main ON/OFF switch

2- Fuse box


3- Take power

4- Voltage-free pedal input

5- 0-10V and 4-20 mA signal input

6- Fan

CONNECTOR (0-10V / 4-20 mA).
0-10V and 4-20 mA CONNECTION




- 1- Positive (+) pin for 0-10V control
- 2- Negative (-) pin for 0-10V control
- 3- not connect.
- 4- Positive (+) pin for 4-20 mA control
- 5- Negative (-) pin for 4-20 mA control



WARNING: Do not use simultaneously the two inputs (0-10 V and 4-20 mA).
It may cause malfunction or breakdown.

3-PIN CONNECTOR (Voltage-free open contact).
ON/OFF CONNECTION



- 1- Pin for connection of voltage-free output (NO). Pump off.
 - 2- Not connected, DO NOT USE.
 - 3- Pin for connection of voltage-free output (NO). Pump off.
- If we physically join pins 1 and 3, the pump will start working.

4.4-FLOW TABLE

Code ▼	Motor rpm	6,4	8,0	9,6	◀ Tube ID mm
1.9749.50	500 ●	220-2700	320-4700	450-6000	Flow ml/min
1.9749.00	330	150-1850	200-2970	250-3300	
● Brushless motor					

Flows calculated with water under normal conditions without outlet back pressure.

5- START UP

Make sure that the mains voltage is between 110 and 230 V.
Connect the power cable to the rear plug and to the main.
Consult the indicative flow table and install the appropriate tube.
See tips in the Tube Change and Head Description section.
Select the desired function.

OPERATING MODES.

Available modes:

- A- Pumping mode.
- B- Ramp mode.
Programming.
Access to the use of the Ramp.
- C- Cyclic dosing mode.
Programming cycles
Access to the use of Cyclic Dosing
- D- Pedal dosing mode.
Programming.
Access to the use of Dosage with pedal.

• Pumping mode.

This mode is used to work with the pump continuously, to transfer liquids.
If the pedal is connected, while it is pressed, the head will work, if it is released, the head will stop.

- 1- Activate rear switch nº1 (photo a) O/1.
- 2- Activate the front switch nº7 (photo a), it will light up in blue and the display will turn on .
- 3- Select the motor speed in % by pressing the decrease or increase keys nº5 or nº6(photo a).
- 4- If it is necessary to change the direction of rotation, press keys no. 1 or no. 2 (photo a).
- 5- Press the No. 4 (photo a) Start / Stop key to start pumping. To stop pumping, use the same key no. 4 (photo a).
If we have the pedal connected, we can start pumping by pressing key no. 4 (figure a) or by keeping the pedal pressed. Working with the pedal, the pump will stop when we stop pressing the pedal.
- 6- The equipment can be programmed so that if there is a power outage, when it is restored, it continues to work at the speed that has been programmed. To program this function and at the speed at which it should resume operation, press key no. 4 (photo a) for 5 seconds. Its pilot light will flash and with keys no. 5 or no. 6 (photo a) we will enter the speed at which we want it to work if there is a power cut and it is restored, it will be memorized by simply pressing key no. 4 (photo a) Start / Stop.

To make loading, purging or cleaning tasks easier, press key no. 3 (photo a) Full. The speed will progressively increase to the maximum and we will see how it increases on the display up to 99. Pressing the Full key again will recover the programmed speed, lowering the display to that speed.

• Ramp mode.

This working mode consists of achieving an increase or decrease in the speed of the spindle, from an initial value to a final value during a set time.

If the first value is less than the final value, there will be an increase; otherwise, if the programmed initial speed is greater than the final one, there will be a decrease.

The time that we can program for the ramp is from 1 minute to 99 minutes.

Programming:

- 1- Activate rear switch nº1 (photo a) O/1.
- 2- Keep keys nº2 and nº3 (photo a) pressed while operating the front switch nº7 (photo a). It will remain illuminated in blue and the yellow led of key no. 2 (photo a) will light intermittently. The displayed value indicates the programmed time for the ramp in minutes. Range: 1 to 99 minutes
- 3- If you want to change, press keys no. 5 or no. 6 (photo a) to set the new time.
- 4- Press key nº4 (photo a) Start / Stop to memorize the chosen value.

- 5- The initial speed of the ramp will then be displayed on the screen.
- 6- If you want to modify it, press keys no. 5 or no. 6 (photo a) to set the new value of the initial speed in %
- 7- Press key nº4 (photo a) to memorize the chosen value.
- 8- The final speed of the ramp will then be displayed on the screen.
- 9- If you want to modify, press the keys nº5 or nº6 (photo a) to set the new value of the final speed in %
- 10- Press key nº4 (photo a) to memorize the chosen value.
- 11- To save all the ramp parameters, press the front switch no. 7 (photo a), the blue light and the equipment display will turn off.

Access to use the Ramp:

From the previous position in which switch no. 7 (photo a) has the blue light off, to access ramp mode and be able to work with the programmed ramp, key no. 2 (photo a) must be held down while the front switch is actuated nº7(photo a). The blue light will turn on, the yellow led of key no. 3 (photo a) will turn on and the display will show the initial speed of the ramp in % and the ramp will start. On the display we will see how the speed increases until reaching the programmed final speed in the time that we have programmed.

All keys will remain inactive during the process.

At the end of the ramp, the spindle will stop. If we want to make another ramp, press any key.

To exit the ramp function, activate the front switch nº7 (photo a). The blue light and the display will go off, and the pump will remain in standby waiting to work in the mode we choose. • **Cyclic dosing mode.**

This work mode is used to be able to work with the pump running for a running time and stopping for a stop time, thus cyclically until we stop the equipment.

Normally it is used for filling a fixed volume of several containers, having a stop time to be able to pass the end of the tube from one container to another.

The time that we can program, both on and off, is from 1 second to 99 seconds.

Before starting the programming, a series of tests must be carried out to verify what volume we want to dose, to know how long the pump must work and at what speed to achieve said volume.

Programming:

- 1- Activate rear switch nº1 (photo a) O/1.
- 2- Keep keys nº1 and nº3 (photo a) pressed while operating the front switch nº7 (photo a). It will remain illuminated in blue and the yellow led of key no. 1 (photo a) will light intermittently. The displayed value indicates the running time in seconds.
Range: 1 to 99 seconds
- 3- If you want to change, press keys no. 5 or no. 6 (photo a) to set the new running time in seconds.
- 4- Press key nº4 (photo a) Start / Stop to memorize the chosen value.
- 5- Then the stopping time will be displayed on the screen.
- 6- If you want to modify it, press keys no. 5 or no. 6 (photo a) to set the new stop time value in seconds.
- 7- Press key nº4 (photo a) to memorize the chosen value.
- 8- Then the motor speed in % will be displayed on the screen.
- 9- If you want to modify it, press keys no. 5 or no. 6 (photo a) to set the new speed value in %
- 10- Press key nº4 (photo a) to memorize the chosen value.
- 11- To store all the parameters of the cyclic dosage, press the front switch no. 7 (photo a), the blue light and the equipment display will turn off.

Access to the cyclic dosage:

From the previous position in which switch no. 7 (photo a) has the blue light off, to access the cyclic dosage mode and be able to work with the programmed cyclic dosage, key no. 1 (photo a) must be held down while it is activated the front switch nº7 (photo a). The blue light will turn on and the yellow led of key no. 2 (photo a) will light up and the display will show the programmed dosing time and the cycle will start. On the display we will see how the time decreases until it reaches zero, then the stop time will appear, and it will go down to zero, like this cyclically until we stop the equipment.

All keys will remain inactive during the process.

To exit the cyclic dosing function, activate the front switch no. 7 (photo a). The blue light and the display will go off, and the pump will remain in standby waiting to work in the mode we choose.

• Pedal dosing mode.

This work mode is used to be able to work with the pump running for a running time when we activate the pedal and stopping when the time reaches zero.

Normally it is used for filling a fixed volume of several containers, having control by means of the pedal of when the head starts up.

The time that we can program is from 1 second to 99 seconds.

Before starting the programming, a series of tests must be carried out to check what volume we want to dose, to know how long the pump must work and at what speed to achieve said volume.

Programming:

- 1- Activate rear switch nº1 (photo a) O/1.
- 2- Keep keys nº1 and nº3 (photo a) pressed while operating the front switch nº7 (photo a). It will remain illuminated in blue and the yellow led of key no. 1 (photo a) will light intermittently. The displayed value indicates the running time in seconds.
Range: 1 to 99 seconds
- 3- If you want to change, press keys no. 5 or no. 6 (photo a) to set the new running time in seconds.
- 4- Press key nº4 (photo a) Start / Stop to memorize the chosen value.
- 5- Then the stopping time will be displayed on the screen.
- 6- If you want to work with the pedal, press keys no. 5 or no. 6 (photo a) to set the stopping time value to 0 seconds.
- 7- Press key nº4 (photo a) to memorize the chosen value.
- 8- Then the motor speed in % will be displayed on the screen.
- 9- If you want to modify it, press keys no. 5 or no. 6 (photo a) to set the new speed value in %
- 10- Press key nº4 (photo a) to memorize the chosen value.
- 11- To store all the parameters of the cyclic dosage, press the front switch no. 7 (photo a), the blue light and the equipment display will turn off.

Access to dosage with pedal.

Connect the pedal to the rear three-pin connector no. 4 (photo b).

From the previous position in which switch no. 7 (photo a) has the blue light off, to access the dosing mode with the pedal and to be able to work with the programmed dosage, key no. 1 (photo a) must be held down while it is activated the front switch nº7 (photo a). The blue light will turn on and the display will show the programmed dosing time.

When we press the pedal, the time will begin to decrease until it reaches zero and the programmed dosage will be carried out. When it reaches zero, the spindle will stop until we press the pedal again.

You must press the pedal, not hold it down.

All keys will remain inactive during the process.

To exit the dosing function with the pedal, activate the front switch no. 7 (photo a). The blue light and the display will go off, and the pump will remain in standby waiting to work in the mode we choose.

NOTE:

While using the rear 0-10V and 4-20mA analog input, the ramp and cyclic dosage modes are not accessible.

The footswitch connection is active during use of the footswitch pump and cycle function.

When using the 0-10V / 4-20mA adjustment, the adjustment keys 5 and 6 are inactive.

6- CHANGE OF TUBES

Press the OFF switch. Extract the tube according to the indications described in the section “253-3r head description” on page 5 of this manual.

When the new tube is installed, it should be centered over the rollers to prevent the rotor from pinching it.

Check that the pump is OFF.

On the sides of the head where the peristaltic tube enters and exits there is a sliding button that acts on the tube fixers. Slide the button to release or clamp the tube. When you release the button it will recover the position by itself.

In general, new tubes can stretch during the first 30 minutes of operation. If this happens, they must be tightened again to avoid unexpected breakage. To detect the elongation of the tube to the head it is useful to conveniently mark the tube with a marker.

Since the friction of the tubes with the rollers increases with the diameter of the tubes, the minimum adjustable speed increases the greater the diameter of the tube or the hardness of the tube.

It is advisable to use a starting speed slightly higher than the minimum, since at any moment it can stop and cause the regulation circuit to overheat, which could be damaged if it remains in this situation for a long time.

When using two heads at the same time, it may be necessary to limit the diameter of the tubes to be used, especially with tubes of high hardness.

The pump supply and discharge tubes can have any wall thickness, but not the tube that is installed in the head, whose wall must be 2,4mm.

The silicone tubes supplied with each pump are medical / food grade according to FDA and USP standards, autoclavable at 120° C, with a peristaltic range of use up to 80° C and medium duration.

Other materials available are:

The most mechanically resistant tubes are PHARMA, TYGON A-60-C ®, TYGON A-60-G ® and medium-duration SILICONE and VITON Ø.

However, the durability also depends to a large extent on the chemical nature of the pumped liquid, the pressure, the existing temperature and naturally the engine revolutions.

Proper choice of tubing ID prevents higher RPM demand from the peristaltic pump motor with small diameter tubing and decreased tubing life.

AVAILABLE MATERIALS:

PHARMA	Autoclavable multiple times. ETO and Gamma sterilizable. Medical-food grade, class VI USP, 21CFR 177.2600 and FDA. Not hemolytic. Excellent resistance to chemicals. ISO 10993. Low permeability and good resistance to abrasion. Long duration. Use temperature, -51°C to 132°C. Beige.
SILICONE	Autoclavable. The most versatile tube. Platinum Cure quality silicone. Average duration. Medical/Food Grade. Excellent biocompatibility. Maximum temperature. 140°C. Translucent.
TYGON A-60-C ®	Autoclavable multiple times. food grade Long duration. Resistant to acids, alkalis, oxidizing agents. Use temperature: -59°C to 135°C. Beige.
TYGON A-60-G ®	Autoclavable Compatible with Ozone, UV light and disinfectants. Great resistance to fatigue and abrasion. Resistant to acids, alkalis and alcohols. Use temperature -59°C to 135°C. Black color.
VITON ®	Autoclavable Suitable for acids and non-acetone solvents. Maximum temperature 300°C. Black color.

2,4 mm wall calibrated tube codes, 1 meter

▼ Tube ID ►	6,4 mm	8,0 mm	9,6 mm
PHARMA	1.8802.64	1.8802.80	1.8802.96
SILICONE	1.8762.64	1.8762.80	1.8762.96
TYGON A-60-G ®	1.8756.64	1.8756.80	1.8756.96
VITON ®	1.8791.64	1.8791.80	1.8791.96

Important: Head tubes should be lightly coated with silicone grease to extend life and ease starting at low rpm.

Silicone grease, 50g for lubricating peristaltic tubes. Code 8.0030.03

7- ORDERING INFORMATION

Code ▼	Motor-rpm	Head	Article
1.9749.50	500	253-3r	Complete Peristaltic Pump
1.9749.00	330	253-3r	Complete Peristaltic Pump
1.0078.44		253-3r	Complementary head
1.9740.01			Foot switch

8- MAINTENANCE-SPARE PARTS

Before any examination or repair of the appliance, it is necessary to disconnect the mains plug. Any initiative must be carried out by qualified personnel to avoid greater evils.

Entrust your device to a technical service authorized by *DINKO Instruments*.



The engine and its block do not require greasing, so they do not have maintenance.

The rotor bearings are self-lubricating, but it is advisable to lightly lubricate them with silicone grease ref. 8.0030.03 or similar from time to time next to the rollers or the head opening lever and its guides, especially if they have been washed.

See Figure.

The head tube must be replaced periodically in a systematic way to avoid the inconvenience of its breaking during full operation of the pump.

Code	-	Description
1.0077.38		Control driver for motor 500 rpm
1.9740.01		Foot switch.
1.0078.44		Head 253-3r
1.0060.30		Main control circuit.
1.0080.01		Motor for head 253-3r, 330 rpm.
1.0077.37		Motor for head 253-3r, 500 rpm
1.8093.21		Power supply 100 W - 24 V.

9- ACCESSORIES

9.1 Balance for flow and dosage calibration.

Reproducibility 0,1 g. 600g capacity. Code 1.9812.02

To measure the quantity dosed in the Calibration process of peristaltic pumps, it is very effective to use a precision balance with digital reading.



If the liquid to be pumped has density "1" there will be no difference between grams and millilitres. Otherwise, calculate the density by weighing a quantity of the liquid with the help of a 25ml test tube, for example, previously taring the test tube on the scale.

Divide the weight indicated on the digital readout of the scale in grams by the millilitres contained in the test tube to obtain the density according to the relationship,

$$D = M / V.$$

There is always the option to Calibrate the pump directly based on weight instead of volume.

Specifications:


- ◆ Single digital reading platter, with highly visible backlit LCD screen.
- ◆ Simple use of great robustness with ABS casing and hermetic anti-humidity membrane keyboard
- ◆ Stainless steel pan, 157x128mm ◆ External auto calibration ◆ Units of measure: grams, pounds and ounces.
- ◆ Continuous tare up to 600 g ◆ Power supply 230V 50/60Hz
- ◆ Non-slip rubber feet
- ◆ Working temperature: from +5°C to +40°C. Maximum use humidity, 85% RH

9.2 Graduated cylinder, 25 ml. Code 1.9808.20

9.3 Silicone grease, 50g. Lubrication of peristaltic tubes. Code 8.0030.03

9.4 Standing support. Code 1.8003.08

Useful as a support for the tube/dispensing tip. Foot: 150 x 70cm. Bar, height 70cm. Sliding support for dosing tip.



CONNECTORS FOR PERISTALTIC TUBES



9.5 Reducing Splice Connectors / Same Ends, polypropylene

For tubes 1,6/3,2mm. ID. Code 1.0080.15

For tubes 3,2/4,8mm ID. Code 1.0080.18

For tubes 4,8/6,4mm ID. Code 1.0080.05

For tubes 6,4/8mm ID. Code 1.0080.14

For tubes 8/12,7mm ID. Code 1.0080.20



9.6 Straight connector for fitting/reducer, polypropylene

Straight connector / reducer Ø 4-5-8 to 7-10-12mm. Code 1.0120.31



9.7 Form Y connectors, polypropylene

Y shape connector, 6mm. either. Code 1.0120.26

Y shape connector, 8mm. either. Code 1.0120.48

Y shape connector, 10mm. either. Code 1.0120.32

Y-shaped connector, 12mm. either. Code 1.0120.33

9.8 Connectors-stainless steel 316 tubes- Connection and dosage

Straight connection tubes 40 mm length

Tube for 0,5 and 0,8mm ID tubes, 25 Units. Code 8.0056.14

Tube for 1,6mm ID tubes, 25 Units. Code 8.0056.06

Tube for 3,2mm ID tubes, 25 Units. Code 8.0056.08

Tube for 4,8mm ID tubes, 25 Units. Code 8.0056.10

Tube for 6,4mm ID tubes, 25 Units. Code 8.0056.12



Dosage 130 mm length with a bevel- Connection and dosage

Dosage tube for 0, 5 and 0, 8 mm ID tubes, 10 Units. Code 8.0056.15

Dosage tube for 1,6 mm ID tubes, 10 Units. Code 8.0056.07

Dosage tube for 3,2 mm ID tubes, 10 Units. Code 8.0056.09

Dosage tube for 4,8 mm ID tubes, 10 Units. Code 8.0056.11

Dosage tube for 6,4 mm ID tubes, 10 Units. Code 8.0056.13



Length 38mm

Tube 0,8mm outer Ø, 10 Units Code 1.0077.23

Tube 0,9mm outer Ø, 10 Units Code 1.0077.26

Clamping flange P. Code 1.0120.01

Clamping flange G. Code 1.0120.12

9.9 Anti-floating sunk Head

For peristaltic tubes of 1,6 /3,2 mm ID. **Code 1.0303.10**

For peristaltic tubes of 4,8 mm.ID **Code 1.0303.11**

For peristaltic tubes of 6,4 mm ID. **Code 1.0303.12**

For peristaltic tubes of 8,0 mm ID. **Code 1.0303.13**

For peristaltic tubes of 9,6 mm ID. **Code 1.0303.14**

For peristaltic tubes of 12,7 mm ID. **Code 1.0303.15**



9.10 Stainless steel tubes for dosing with backflow valve

For tubes size 3,2/4,8 mm. ID. Tip 4 mm. ID. 1mm.TW. **Code 1.0302.10**

For tubes size 4,8/6,4 mm. ID. Tip 6 mm. ID. 1 mm.TW. **Code 1.0302.11**

For tubes size 6,4/8 mm. ID. Tip 8 mm. ID. 1 mm. TW. **Code 1.0302.12**

For tubes size 8/9,6 mm. ID. Tip 10 mm. ID. 1 mm. TW. **Code 1.0302.13**



10- CHANGE OF FUSES

The fuse box is part of the power base located at the rear of the pump. See Figure.



Main switch

Fuse box

Power base

Pry with a screwdriver between the central part of the fuse holder box and the upper part of the power supply base to remove the fuse holder box. The box remains attached without being fully extracted. There are two fuses.

Press the box in to restore its original position.

Remember that you have already used the spare fuse.

11- WARRANTY

11.1 DURATION:

The guarantee is established for a period of 1 year from the date of commissioning of the device, provided that the guaranteed card is returned to us within 8 days of said commissioning.

Without this condition the guarantee will not be valid.

11.2 SCOPE OF WARRANTY:

The guarantee is given against manufacturing and material defects for an average work week of 40 hours.

The guarantee is reduced proportionally to the increase in working hours.

Repairs will be made in our factory.

Otherwise, the guarantee will only include the replacement of the defective elements.

DINKO will not be responsible for the transportation costs, nor will it assume responsibility for the consequences caused by the immobilization of the device.

The parts replaced free of charge remain our property, reserving the right to request their return, free of postage to our address.

Repairs or replacement of parts during the warranty period do not extend the initial warranty.

Our responsibility is limited to the attached guarantee and not to possible accidents to persons or other things.

Any alteration of the device by the user voids the guarantee.

12- "CE" DECLARATION OF CONFORMITY

DINTER S.A - *DINKO Instruments* c/ Encarnació, 123-125 / 08024- Barcelona

Declares that the items mentioned in the attached list, to which this declaration refers, comply with the essential safety requirements of the applicable European Directive:

- Low Voltage Directive Directive D2006/95/CEE of December 12, 2006
- Essential requirements of Annex I of the Machinery Directive 2006/42/CEE of May 17 from 2006
- Electromagnetic Compatibility Directive 2004/108/CEE of December 15, 2004
- Safety for electrical measurement, control and laboratory devices. Requirements relating to the EMF. IN 61326
- Safety rules for electrical measurement, control and laboratory devices. Part I. General prescriptions EN 61010-1

However, the user must observe the assembly and connection instructions indicated in the technical instructions catalogues.

Name	Joan A. Bravo	Josep X. Sensada
Position:	Technical Director	Responsible for Quality

Signature



Model: Peristaltic Pumps D-25Vplus 1.9749.00 and 1.9749.50

OTHER DINKO APPARATUS

- Blenders-Homogenizers
- Colorimeters
- Conductivity Meters
- Dosing Pumps
- Extractor for meat analysis
- Heating Plates
- Infrared Stoves
- Kits for water analysis
- Magnetic Stirrers
- Metallic block heaters
- Microscopes
- Nephelometers
- Orbital Shakers
- Oximeters
- Peristaltic Pumps
- pH-meters
- Photometers
- Respirometers
- Rod Stirrers
- Rotary Stirrers
- Sand Baths
- Spectrophotometers
- Temperature Controllers
- Timers / *Timers*
- Trichinoscope - TriquiVisor
- Turbidity Meters
- Turn dishes
- Vacuum Pumps



DINTER, S.A c/ Encarnació, 123-125. Tel. +34 93 284 69 62 - 08024-Barcelona

dinter@dinko.es www.dinko.es