

INFRA - RED OVEN

Model D - 60

The DINKO Infra-red Oven is constructed for a quick heating of samples or containers from all the possible angles to a closely controlled desired temperature under observation

Oven has louvers for exterior cooling and are finished in acid resistant, brown and beige baked enamel.

The heat reflecting oven inside walls mirror finished stainless steel reflect the heating rays onto the object

The upward sliding, heat resistant, glass door, provides a observation window

Temperatures of 200°C are obtained in about ten minutes, for what it is not necessary that the ovens remains always connected



The control panel contains general switch, electronic temperature controller with digital indication and lamp synchronised

Specifications

Digital Temperature Control: maxim 200 °C.

Accuracy: $\pm 1^{\circ}$ C.

Size: overall 25 x 26 x 40cm. Chamber: 18 x 18 x 11cm.

Weight: 10Kg.

Operate on: 230V 50/60Hz. Power: 350W.

Code	Article
1.9766.00	Infra-Red Oven D-60, digital temperature controller
6.9763.00	Disposable pans for solids, in aluminium, 100pcs
6.9764.00	Disposable pans for liquids, in aluminium, 25pcs

Distributor: