

## COOLER BATH

### COOLER (Peltier) – HEATER THERMOSTATED UNIT WITH RECIRCULATION Model D-40

- ◆ Microprocessor control
- ◆ Stainless steel made
- ◆ Range of temperature control: 2,5 to +60°C
- ◆ Digital indicator
- ◆ Peristaltic pump for recirculation fluid: 4 to 20l/hour to thermostat outside devices
- ◆ Peltier cooling system
- ◆ Removable stainless steel reservoir

Supplied with cover reservoir and tubes



Parameters	Features	Parameters	Features
Temperature range	2,5°C to +60°C	Power	300W (230V) - 600W (110V)
Fluctuation of temperature	± 0.1°C	Heating unit	Cartridge
Heating power	100W	Recirculation control	from 4 to 20L/hour
Cooling capacity	180W	Reservoir capacity	1L /10x15cm, Ø x height
Coolant	H <sub>2</sub> O	Overall size	27x34x29cm width x deep x height
Heating time	<30 min from 20°C to 60°C	Operation/ Weight	110-230V 50/60Hz / 9.5Kg
Cooling time	<100 min from 20°C to 5°C	Code	1.8200.00*

Distributor: