

BLENDER. Model CB15

High capacity, 4 litres, power and speed. 3.75 HP motor
 It crushes, emulsifies, mixes, pulverizes, homogenizes, or disintegrates in seconds
 Stainless steel container with blades, independent of the base motor
 Adapting AD1 to use containers of 1L and smaller
 Three speeds or continuous speed control, models CB15VE and CB15VEX
 Operation 230V 50/60Hz. 2250W

High speed, 20.000rpm free, 15.500rpm with 2L of water, 14.000rpm with 4L
 Medium speed 18.300rpm free, 14.400rpm with 2L of water, 13.200rpm with 4L
 Low speed 15.500rpm free, 12.500rpm with 2L of water, 11.500rpm with 4L



Code	Articles
10.9797.00	Blender CB15E, 4L, with lid, three speeds
10.9798.05	Blender CB15TE, 4L, with lid, three speeds and 3 minutes digital timer
10.9798.00	Blender CB15VE, 4L, with lid, continuous speed analogic control
10.9798.06	Blender CB15VEX, 4L, with lid, continuous speed digital control
10.9799.00	Adapter AD1 for use 1L, 1.25 L container and mini MC containers

Blender. Model 8010

High power, 400W
 Two speeds, 22000rpm and 18500rpm
 One litre stainless steel or 1,25 litre glass container with vinyl lid
 One litre and semimicro containers with screw cover all stainless steel.
 Blades integrated in the containers
 Independent container of the base - motor
 Mini containers stainless steel in option with silicone lid
 Timer 1 to 180 seconds. Operation 240V 50/60Hz



Base Blender only



Mini containers MC1, MC2 and MC3



360 ml



One litre with screw lid



Hermetic container



1L st. steel and 1.25L glass



Assembly, axis / blades



Set of tools

Code	Articles
10.9794.00	Blender 8010ES. One litre stainless steel container. Two speeds. 2 pieces vinyl lid
10.9795.00	Blender 8010EG. 1.25L glass container. Two speeds. 2 pieces vinyl lid
10.9794.17	Blender only base
10.9794.01	One litre stainless steel container 2-piece vinyl lid, spare
10.9795.01	1,25L glass container 2-piece vinyl cover, spare
10.9807.00	One litre container with screw cover and gaskets, spare. Working from 75 to 500 ml. All stainless steel
10.9794.16	One litre container, hermetic cover and sample valve. All stainless steel
10.9800.00	Mini container MC1, stainless steel, 12-37 ml. Silicone lid
10.9801.00	Mini container MC2, stainless steel, 37-110 ml. Silicone lid
10.9802.00	Mini container MC3, stainless steel, 50-250 ml. Silicone lid
10.9806.00	Semi-micro container 360 ml. Working from 25 to 250 ml. Screw cover and gaskets. All stainless steel
10.9794.10	Blades for 8010
10.9800.01	Blades for MC1
10.9801.01	Blades for MC2
10.9802.01	Blades for MC3
10.9794.02	Assembly, axis / blades / gasket PTFE
10.9794.18	Set of tools to remove blades

Distributor: