

## MAGNETIC STIRRER 4 PLACES

### Model D-66AM, with metallic block

Magnetic stirrer incubator for 4 places of flasks, tubes, etc. of 40mm. diameter.

Available aluminium adapters for smaller containers of 40mm diameter.

Deep holes, 60mm

Temperature control from room to 50°C

Independent magnetic stirring for each place



The front panel include the temperature display of the PID touch Controller to temperature set point, 4 light stirring switches for each position of the metal block and general pilot light.

Suitable for incubation tasks that require constant agitation to process acceleration

The rear panel include the ON-OFF switch and the main connector with fuse.

Supplied with 4 magnetic bars 2x7mm or 5x12mm

### Heather block Model D-66

The D-66 incubator is like to model D-66AM but without magnetic stirrer

#### Specifications:

Temperature control: room to 50°C.

Precision:  $\pm 0,1^{\circ}\text{C}$

Speed: 3000rpm

Size: 204 x 285 x 130mm

Weight: Model D-66AM, 4,8Kg. Model D-66, 3,8Kg

Operation: 230VAC 50/60 Hz. 200W

Code	Model	Article
1.8093.00	D-66AM	Magnetic stirrer 4 places with metallic heather block
1.8093.20	D-66	Heather block incubator 4 places
1.8093.15		Adapter of aluminum (to indicate inside diameter)
1.8049.00		Magnetic bar 2x7mm
1.8050.00		Magnetic bar 5x12mm

Distributor: