

PERISTALTIC PUMPS, MULTIHEADS, VARIABLE FLOWS Models D25V2i, D25V4i and D25VT2i

- ◆ Numeric speed control that can be set 0 - 99% in 1% increments with 100% repeat accuracy by pump head.
- ◆ Remote control by foot switch or external relay
- ◆ Pump heads 50-2r of 2 rollers or 4 rollers for CFV-4r
- ◆ Switch ON/OFF/Reverse flow by pump head
- ◆ Digital cyclic timer allows pumping schedule and shutdown times of 0,01 seconds at 99,90 hours which makes it possible to obtain any volume in each cycle of the timer
- ◆ Speed control by optional signal 0-10 volts. Consult this option.
- ◆ Pumps accept several tube sizes. Tubing is easily removed for changeover, or autoclaving.
- ◆ Pumps with multiple heads adapted to the needs of the user are manufactured.



D-25V2i




D-25V4i

Ordering information

Code ▼	Model	Pump heads	Cyclic timer	Flow ml/min	Motor rpm
1.9733.05	D25V2i	2x50-2r	No	0,1-100	30
1.9738.00	D25V2i	2x50-2r	No	0,4-300	80
1.9738.01	D25V2i	2x50-2r	No	1-700	240
1.9734.05	D25V4i	4x50-2r	No	0,1-100	30
1.9739.00	D25V4i	4x50-2r	No	0,4-300	80
1.9739.01	D25V4i	4x50-2r	No	1-700	240
1.9733.11	D25VT2i	2xCFV-4r	Yes	0,02-25	15
1.9733.12	D25VT2i	2xCFV-4r	Yes	0,03-30	30
1.9733.04	D25VT2i	2xCFV-4r	Yes	0,08-85	80
1.9738.17	D25VT2i	2xCFV-4r	Yes	0,25-250	240

PUMP HEAD CFV-4r



- 1- Mobile head cover for tube change
- 2- Fixed lower cover
- 3- Tube retainer
- 4- Peristaltic tube
- 5- Release lever for head change

Size: with 2 pump heads, 32 x 26 x 14cm. Weight: 5 / With 4 pump heads, 26 x 34 x 27cm. Weight: 6Kg
Operation: 100... 240V 50/60Hz

Distributor: