

MULTI-HEAD PERISTALTIC PUMPS - VARIABLE FLOW

Models D25V2i, D25V4i and D25VT2i or D25VT4i with cyclic timer.

- ◆ Independent numerical speed control for each head with 1% repeatability.
- ◆ Remote control by foot pedal or external relay.
- ◆ 4 easy-loading pump heads available with 3 or 4 roller rotors and gear motors of different speeds
- ◆ Individual ON/OFF switch for each head with flow inverter.
- ◆ Digital cyclic timer on T models to program repetitive pumping and stopping times. 0.01 seconds to 99.90 hours, making it possible to obtain any volume in each cycle.
- ◆ Optional 0-10 volt analogue signal speed control. Consult this option.

Pumps with multiple pump heads are manufactured to suit the user's needs.

PUMP HEADS



CFG1-3r head, 1 channel



CFG2-3r head, 2 channels



CFV-4r head



CF-4r head

CFG heads are opened by turning the front knob, while CFV and CF heads are opened by pulling up on the top, which makes it easier to remove or install the peristaltic tubes.

Pump Codes. 100... 240V 50/60Hz. Indicative flow rates per channel

Code ▼	Model	Pump Head	Timer	Flow ml/min	Engine rpm
1.9733.14	D25V2i	2 x CFG1-3r	No	5.0 to 565	300
1.9733.15	D25V2i	2 x CFG2-3r		0.5 to 175	300
1.9733.40	D25V2i	2 x CF-4r		2.0 to 1100	350
1.9733.13	D25V2i	2 x CFV-4r		0.2 to 30	30
				0.5 to 80	80
				2.0 to 220	240
1.9733.16	D25VT2i	2 x CFG1-3r	Yes	5.0 to 565	300
1.9733.17	D25VT2i	2 x CFG2-3r		0.5 to 175	300
1.9733.41	D25VT2i	2 x CF-4r		2.0 to 1100	350
1.9733.04	D25VT2i	2 x CFV-4r		0.2 to 30	30
				0.5 to 80	80
				2.0 to 220	240
1.9734.01	D25V4i	4 x CFG1-3r	No	5.0 to 565	300
1.9734.02	D25V4i	4 x CFG2-3r		0.5 to 175	300
1.9734.08	D25V4i	4 x CFV-4r		0.2 to 30	30
				0.5 to 80	80
				2.0 to 220	240
1.9734.11	D25VT4i	4 x CFG1-3r	Yes	5.0 to 565	300
1.9734.12	D25VT4i	4 x CFG2-3r		0.5 to 175	300
1.9734.10	D25VT4i	4 x CFV-4r		0.2 to 30	30
				0.5 to 80	80
				2.0 to 220	240

Sizes: With 2 heads, 32 x 26 x 14 cm. Weight: 5 kg. With 4 heads, 26 x 34 x 27 cm. Weight: 7 kg .

Distributor: