

ROTATORS

ROTATORS FOR TUBES or CAPSULES

Model D- 01

Rotator provided of disks with clamps for tubes or tray for capsules
Angle of rotation can be continuously set between 0 and 90°

Different disks available allow the use of tubes of different sizes, from Eppendorf tubes to tubes of 24 mm in diameter
Ask other options

Standard trays offer for capsules or bottles of 28 or 40mm diameter
Ask other options

Comes with line cord but without disc or tray

Speed control: 0-15rpm. 230V 50/60Hz.



Code	Article
1.8007.00	Rotator stirrer D-01
1.8016.00	Tray for 18 capsules to 28mm Ø
1.8017.00	Tray for 8 capsules to 40mm Ø
1.8016.01	Tray mix (to indicate)
1.8008.01	Disk for 16 Eppendorf tubes
1.8008.00	Disk for 16 tubes of Ø 12/14mm
1.8009.00	Disk for 16 tubes of Ø 15/18mm
1.8010.00	Disk for 16 tubes of Ø 20/24mm

Distributor: