

Model TriquiVisor Automat XY- HD

The *Dinko* TrichiVisor Automat XY- HD, it is a new instrument for a sequential, precise and quick exam of tissue animals specially conceived for the detection of trichinas

This instrument is according of the official method for the identification of trichinas published in the Official Newspaper of European Union August 11, 2015, Implementing Regulations (EU) 2015/1375 of the Commission on August, 10, 2015.

Features:

- Panoramic view finder of high resolution HD and contrast that allows to see a whole grid of the collection cuvette shown with 40X magnification
- Precise and quick focusing system with helical ring
- White LED lighting crown with evenly distributed power equivalent to 100 watts
- Filter diffuser with adjustable light polarizer which prevents annoying reflections for viewing or loss of contrast
- Scale continues selectable from 35X to 160X magnification.
- Moving motorized horizontal and vertical viewing independently with type switches pushbutton start / stop
- Selectors for independent horizontal and vertical displacements with three positions: right -anchor - left for horizontal and up - anchor - down to the vertical that guarantee systematic and sequential observation grids of the sample cuvette
- Sample cuvette - holder board for stability, with grids
- Protocol for centering the test image using a 90X
- It meets the essential safety requirements of the European Directive applicable (CE)

TriquiVisor Automat XY - HD: W46 x D40 x H28cm. Code 1.9810.00, screen LED 24" and cuvette included
TriquiVisor Automat XY - HD: W46 x D40 x H28cm. For PC. Code 1.9810.30, cuvette included, without PC
Weight: 18Kg
Operation: 230V 50/60 Hz.

Screen: LED high resolution 24". Code 1.9810.23
Acrylic cuvette 180 x 40mm. Code 1.9810.22
PC with 24" LED screen. Code 8.9749.00

Distributor:

