

TURN DISHES

MODEL D-10

Apparatus provided of a round platform that can be turn slowly to 6 r. p. s.

On front panel can be find the speed control button, pilot lamp and ON / OFF switch.



The platform have round marks for to place Petri Dishes of 55, 90, 100 and 140mm diameter.

The rotation to make easy a uniform micro-organism scatter on Petri dishes.

On rear panel the Turn Dishes are provided of foot switch connector. The foot switch is supplied as accessory and to make possible the start / stop function for hand free.

Speed: 0 - 360 rpm.

Overall size: 160 x 200 x 120mm

Platform size: 150mm diameter.

Weight: 1,5kg.

Operate: 230V 50/60Hz. 0,3Amp.

Code	Articles
1.8018.00	Turn Dishes D-10
1.9740.00	Foot switch

Distributor: